

MBD4 Antibody

PACO10407

Description

This MBD4 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO10407
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti human PolyClonal Antibody
Synonyms:	3 N(4) ethenocytosine glycosylase antibody, G/5 fluorouracil mismatch glycosylase with biphasic kinetics antibody, G/T mismatch glycosylase antibody, G/U mismatch glycosylase antibody, MBD 4 antibody, MBD4 antibody, MBD4_HUMAN antibody, MED 1 antibody, MED1 antibody, Methyl CpG binding domain protein 4 antibody, Methyl CpG binding endonuclease 1 antibody, Methyl CpG binding protein MBD4 antibody, Methyl-CpG-binding domain protein 4 antibody, Methyl-CpG-binding endonuclease 1 antibody, Methyl-CpG-binding protein MBD4 antibody, Mismatch specific DNA N glycosylase antibody, Mismatch-specific DNA ... Read more
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse

Antibody Data

Isotype:	IgG
Uniprot:	O95243
Host Species:	Rabbit
Purification:	Antigen Affinity purified
Immunogen:	Human MBD4
Immunogen Species:	Homo sapiens (Human)

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Buffer: PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Form: -

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-