

MSL3 Antibody

PACO10624

Description

This MSL3 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO10624
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti human PolyClonal Antibody
Synonyms:	MSL3 antibody, AU018931 antibody, DKFZp586J1822 antibody, Drosophila MSL3 like 1 antibody, Male specific lethal 3 (Drosophila) like 1 antibody, Male specific lethal 3 homolog (Drosophila) antibody, Male specific lethal 3 homolog 1 (Drosophila) antibody, Male specific lethal 3 homolog 1 antibody, Male specific lethal 3 homolog antibody, Male specific lethal 3 like 1 (Drosophila) antibody, Male specific lethal 3 like 1 antibody, Male specific lethal 3 protein like 1 antibody, Male-specific lethal 3 homolog antibody, Male-specific lethal-3 homolog 1 antibody, Male-specific lethal-3 protein-like 1...Read more
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q8N5Y2
Host Species:	Rabbit
Purification:	Antigen Affinity Purified
Immunogen:	Human MSL3
Immunogen Species:	Homo sapiens (Human)

Buffer: PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Form: -

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-