

NMD3 Antibody

PACO10872

Description

This NMD3 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO10872
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti human PolyClonal Antibody
Synonyms:	60S ribosomal export protein NMD3 antibody, CGI 07 antibody, FLJ21053 antibody, hNMD3 antibody, NMD3 antibody, NMD3 homolog antibody, NMD3 homolog (S. cerevisiae) antibody, NMD3 ribosome export adaptor antibody, NMD3, S. cerevisiae, homolog of antibody, NMD3_HUMAN antibody
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q96D46
Host Species:	Rabbit
Purification:	Antigen Affinity purified
Immunogen:	Human NMD3
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form:	-

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-