

PCDH11Y Antibody

PACO11094

Description

This PCDH11Y Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO11094
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti Human PolyClonal Antibody
Synonyms:	PCDH11Y antibody, PCDH11 antibody, PCDH22 antibody, PCDHY antibody, Protocadherin-11 Y-linked antibody, Protocadherin-11 antibody, Protocadherin on the Y chromosome antibody, PCDH-Y antibody, Protocadherin prostate cancer antibody, Protocadherin-PC antibody, Protocadherin-22 antibody
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q9BZA8
Host Species:	Rabbit
Purification:	Antigen Affinity purified
Immunogen:	Human PCDH11Y
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form:	-

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-