

PLEKHF1 Antibody

PACO11299

Description

This PLEKHF1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO11299
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti human PolyClonal Antibody
Synonyms:	PLEKHF1 antibody, APPD antibody, LAPF antibody, ZFYVE15 antibody, Pleckstrin homology domain-containing family F member 1 antibody, PH domain-containing family F member 1 antibody, Lysosome-associated apoptosis-inducing protein containing PH and FYVE domains antibody, Apoptosis-inducing protein antibody, PH and FYVE domain-containing protein 1 antibody, Phafin-1 antibody, Zinc finger FYVE domain-containing protein 15 antibody
Clone:	-
Applications:	ELISA WB IF
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q96S99
Host Species:	Rabbit
Purification:	Antigen Affinity purified
Immunogen:	Human PLEKHF1
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Form: -

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-