

PPP1R2P9 Antibody

PACO11424

Description

This PPP1R2P9 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO11424
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti human PolyClonal Antibody
Synonyms:	PPP1R2C antibody, PPP1R2P9 antibody, Protein phosphatase inhibitor 2 family member C antibody, PPP1R2 family member B antibody, Protein phosphatase 1 antibody, regulatory subunit 2 pseudogene 9 antibody, Type-1 protein phosphatase inhibitor 4 antibody, I-4 antibody
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	O14990
Host Species:	Rabbit
Purification:	Antigen Affinity Purified
Immunogen:	Human PPP1R2P9
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form:	-

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-