

PSAPL1 Antibody

PACO11515

Description

This PSAPL1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO11515
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti Human PolyClonal Antibody
Synonyms:	Proactivator polypeptide like 1 antibody, Prosaposin like 1 antibody, Prosaposin like protein 1 antibody, PSAPL1 antibody, SAPL1_HUMAN antibody, Saposin A like antibody, Saposin B like antibody, Saposin B Val like antibody, Saposin C like antibody, Saposin D like antibody, Saposin D-like antibody
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q6NUJ1
Host Species:	Rabbit
Purification:	Antigen Affinity purified
Immunogen:	Human PSAPL1
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form:	-

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-