

RNPS1 Antibody

PACO11882

Description

This RNPS1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO11882
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti human PolyClonal Antibody
Synonyms:	E5.1 antibody, MGC117332 antibody, RNA binding protein S1 serine rich domain antibody, RNA binding protein with serine rich domain 1 antibody, RNA-binding protein with serine-rich domain 1 antibody, RNPS 1 antibody, rnps1 antibody, RNPS1_HUMAN antibody, SR protein antibody, SR related protein LDC2 antibody, SR-related protein LDC2 antibody
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q15287
Host Species:	Rabbit
Purification:	Antigen Affinity Purified
Immunogen:	Human RNPS1
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form:	-

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-