

SALL2 Antibody

PACO12018

Description

This SALL2 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

| | |
|----------------------|--|
| SKU: | PACO12018 |
| Contents: | 100µl Bradford Reagent: 1 vial (2ml) |
| Category: | Purified Rabbit Anti human PolyClonal Antibody |
| Synonyms: | SALL2 antibody, KIAA0360 antibody, SAL2 antibody, ZNF795 antibody, Sal-like protein 2 antibody, Zinc finger protein 795 antibody, Zinc finger protein SALL2 antibody, Zinc finger protein Spalt-2 antibody, Sal-2 antibody, hSal2 antibody |
| Clone: | - |
| Applications: | ELISA WB |
| Conjugation: | Non-conjugated |
| Reactivity: | Human, Mouse, Rat |

Antibody Data

| | |
|---------------------------|--|
| Isotype: | IgG |
| Uniprot: | Q9Y467 |
| Host Species: | Rabbit |
| Purification: | Antigen Affinity Purified |
| Immunogen: | Human SALL2 |
| Immunogen Species: | Homo sapiens (Human) |
| Buffer: | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Form: | - |

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

| Recommended Dilutions: | Application | Recommended Dilution |
|------------------------|-------------|----------------------|
| | - | - |

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

| Image | Description |
|-------|-------------|
| - | - |