

TUBA8 Antibody

PACO12984

Description

This TUBA8 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO12984
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Polyclonal Antibody
Synonyms:	Alpha tubulin 8 antibody, Alpha-tubulin 8 antibody, TBA8_HUMAN antibody, Tuba8 antibody, TUBAL2 antibody, Tubulin alpha 8 chain antibody, Tubulin alpha chain like 2 antibody, Tubulin alpha chain-like 2 antibody, Tubulin alpha-8 chain antibody
Clone:	-
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q9NY65
Host Species:	-
Purification:	Antigen Affinity purified
Immunogen:	Human TUBA8
Immunogen Species:	Homo sapiens (Human)
Buffer:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-