

TXNRD3 Antibody

PACO13006

Description

This TXNRD3 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO13006
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti Human PolyClonal Antibody
Synonyms:	TGR antibody, Thioredoxin and glutathione reductase antibody, Thioredoxin glutathione reductase antibody, Thioredoxin reductase 2 antibody, Thioredoxin reductase 3 antibody, Thioredoxin reductase TR2 antibody, TR2 antibody, TRXR3 antibody, TRXR3_HUMAN antibody, TXNRD3 antibody
Clone:	-
Applications:	ELISA WB IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q86VQ6
Host Species:	Rabbit
Purification:	Antigen Affinity purified
Immunogen:	Human TXNRD3
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form:	-

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-