

UNKL Antibody

PACO13097

Description

This UNKL Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO13097
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti human PolyClonal Antibody
Synonyms:	C16orf28 antibody, FLJ12623 antibody, FLJ23360 antibody, KIAA0734 antibody, MGC5179 antibody, Putative E3 ubiquitin-protein ligase UNKL antibody, RING finger protein unkempt-like antibody, Unkempt homolog (Drosophila) like antibody, Unkempt like antibody, Unkempt like protein antibody, Unkl antibody, UNKL_HUMAN antibody, ZC3H5L antibody, ZC3HDC5L antibody, Zinc finger CCCH domain-containing protein 5-like antibody
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q9H9P5
Host Species:	Rabbit
Purification:	Antigen Affinity purified
Immunogen:	Human UNKL
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: -

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-