

XYLT1 Antibody

PACO13270

Description

This XYLT1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO13270
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti Human PolyClonal Antibody
Synonyms:	Beta D xylosyltransferase 1 antibody, O xylosyltransferase 1 antibody, Peptide O-xylosyltransferase 1 antibody, PXYLT1 antibody, UDP-D-xylose:proteoglycan core protein antibody, XT-I antibody, XT1 antibody, XTI antibody, Xylosyltransferase 1 antibody, Xylosyltransferase I antibody, xylosyltransferase iota antibody, Xylt-I antibody, XYLT1 antibody, XYLT1_HUMAN antibody
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q86Y38
Host Species:	Rabbit
Purification:	Antigen Affinity purified
Immunogen:	Human XYLT1
Immunogen Species:	Homo sapiens (Human)
Buffer:	0.1M NaHCO ₃ , 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form:	-

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-