

PLA2G12B Antibody

PACO13566

Description

This PLA2G12B Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO13566
Contents:	100µl Bradford Reagent: 1 vial (2ml)
Category:	Purified Rabbit Anti Human PolyClonal Antibody
Synonyms:	PLA2G12B antibody, PLA2G13 antibody, FKSG71 antibody, Group XIIIB secretory phospholipase A2-like protein antibody, Group XIII secretory phospholipase A2-like protein antibody, GXIII sPLA2-like antibody, sPLA2-GXIIB antibody, GXIIB antibody
Clone:	-
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse

Antibody Data

Isotype:	IgG
Uniprot:	Q9BX93
Host Species:	-
Purification:	Antigen Affinity purified
Immunogen:	Human PLA2G12B-Specific
Immunogen Species:	Homo sapiens (Human)
Buffer:	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form:	-

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-