

PGBD4 Antibody

PACO15205

Description

This PGBD4 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO15205
Contents:	50µl Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	FLJ32638 antibody, FLJ37497 antibody, PGBD4 antibody, PGBD4_HUMAN antibody, PiggyBac transposable element-derived protein 4 antibody
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q96DM1
Host Species:	Rabbit
Purification:	Antigen affinity purification
Immunogen:	Fusion protein of Human PGBD4
Immunogen Species:	Homo sapiens (Human)
Buffer:	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	ELISA	1:1000-1:2000
	WB	1:200-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane 1-2: A549 cells, A375 cells, Primary antibody: PACO15205 at dilution 1/250, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 5 seconds