

DEFA6 Antibody

PACO15816

Description

This DEFA6 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO15816
Contents:	50µl Bradford Reagent: 1 vial (2ml)
Category:	Polyclonal Antibody
Synonyms:	alpha 6 antibody, DEF6 antibody, DEF6_HUMAN antibody, DEFA 6 antibody, DEFA6 antibody, Defensin 6 antibody, Defensin alpha 6 antibody, Defensin alpha 6 Paneth cell specific antibody, Defensin antibody, Defensin-6 antibody, HD 6 antibody, HD6 antibody
Clone:	-
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q01524
Host Species:	Rabbit
Purification:	Antigen affinity purification
Immunogen:	Fusion protein of Human DEFA6
Immunogen Species:	Homo sapiens (Human)
Buffer:	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.


Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	ELISA	1:2000-1:5000
	IHC	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using PACO15816 at dilution 1/40, on the right is treated with fusion protein. (Original magnification: ×200)