

SLC2A12 Antibody

PACO18006

Description

This SLC2A12 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO18006
Contents:	50µl Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Glucose transporter type 12 antibody, GLUT 12 antibody, GLUT-12 antibody, GLUT12 antibody, GLUT8 antibody, GTR12_HUMAN antibody, SLC2A12 antibody, Solute carrier family 2 (facilitated glucose transporter) member 12 antibody, Solute carrier family 2 facilitated glucose transporter member 12 antibody, Solute carrier family 2 member 12 antibody, Solute carrier family 2, facilitated glucose transporter member 12 antibody
Clone:	-
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q8TD20
Host Species:	Rabbit
Purification:	Antigen affinity purification
Immunogen:	Synthetic peptide of Human SLC2A12
Immunogen Species:	Homo sapiens (Human)
Buffer:	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	ELISA	1:1000-1:2000
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

The image is immunohistochemistry of paraffin-embedded Human tonsil tissue using PACO18006 at dilution 1/30. (Original magnification: ×200)