

SYT2 Antibody

PACO18437

Description

This SYT2 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO18437
Contents:	50µl Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Synaptotagmin 2 antibody, Synaptotagmin II antibody, Synaptotagmin-2 antibody, Synaptotagmin2 antibody, SynaptotagminII antibody, SYT 2 antibody, Syt II antibody, Syt2 antibody, SYT2_HUMAN antibody, SytII antibody
Clone:	-
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q8N9I0
Host Species:	Rabbit
Purification:	Antigen affinity purification
Immunogen:	Synthetic peptide of Human SYT2
Immunogen Species:	Homo sapiens (Human)
Buffer:	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	ELISA	1:2000-1:5000
	WB	1:500-1:2000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Gel: 10%SDS-PAGE, Lysate: 50 µg, Lane 1-2: Human fetal brain tissue, Human lung tissue, Primary antibody: PACO18437 at dilution 1/700, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 40 seconds