

ZBTB2 Antibody

PACO20943

Description

This ZBTB2 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO20943
Contents:	50µl Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	bA351K16.2 antibody, Gm1103 antibody, KIAA1483 antibody, ZBTB2 antibody, ZBTB2_HUMAN antibody, Zinc finger and BTB domain containing 2 antibody, Zinc finger and BTB domain-containing protein 2 antibody, ZNF437 antibody
Clone:	-
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q8N680
Host Species:	Rabbit
Purification:	Antigen affinity purification
Immunogen:	Synthetic peptide of Human ZBTB2
Immunogen Species:	Homo sapiens (Human)
Buffer:	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	ELISA	1:2000-1:10000
	IHC	1:30-1:150

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using PACO20943 at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)

The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using PACO20943 at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)