

ALB Antibody

PACO24656

Description

This ALB Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO24656
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	ALB antibody, Serum albumin antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Rabbit

Antibody Data

Isotype:	IgG
Uniprot:	P49065
Host Species:	Guinea pig
Purification:	>95%, Protein G purified
Immunogen:	Serum Albumin purified from Rabbit serum (Native Protein)
Immunogen Species:	Rabbit
Buffer:	Preservative: 0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2
Form:	Liquid

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-