

Abcb4 Antibody

PACO24684

Description

This Abcb4 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO24684
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Abcb4 antibody, Mdr2 antibody, Pgy-2 antibody, Pgy2 antibody, Phosphatidylcholine translocator ABCB4 antibody, EC 7.6.2.1 antibody, ATP-binding cassette sub-family B member 4 antibody, Multidrug resistance protein 2 antibody, Multidrug resistance protein 3 antibody, P-glycoprotein 2 antibody, P-glycoprotein 3 antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Mouse

Antibody Data

Isotype:	IgG
Uniprot:	P21440
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Mouse Phosphatidylcholine translocator ABCB4 protein (352-708AA)
Immunogen Species:	Mus musculus (Mouse)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-