

LGB Antibody

PACO25284

Description

This LGB Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO25284
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	Polyclonal Antibody
Synonyms:	Beta lactoglobulin B antibody, Beta-lactoglobulin antibody, Beta-LG antibody, BLG antibody, LACB_BOVIN antibody, LGB antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Bovine

Antibody Data

Isotype:	IgG
Uniprot:	P02754
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Bovine milk Beta-lactoglobulin (Native Protein)
Immunogen Species:	Bovine
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western blot All lanes: Bovine milk Beta-lactoglobulin antibody at 1:500 Lane 1: Bovine milk Beta-lactoglobulin at 0.1mg/ml Lane 2: Bovine milk Beta-lactoglobulin at 10µg/ml Secondary Goat polyclonal to rabbit IgG at 1/15000 dilution Predicted band size: 18 kDa Observed band size: 18 kDa