

C5 Antibody

PACO25384

Description

This C5 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO25384
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	C5 antibody, Complement C5a anaphylatoxin antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Pig

Antibody Data

Isotype:	IgG
Uniprot:	P01032
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Pig Complement C5a anaphylatoxin protein (1-74AA)
Immunogen Species:	Pig
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:10-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human kidney tissue using PACO25384 at dilution of 1:10