

CRYZL1 Antibody

PACO25961

Description

This CRYZL1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO25961
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	4P11 antibody, Crystallin zeta (quinone reductase) like 1 antibody, CRYZL1 antibody, Protein 4P11 antibody, QOH 1 antibody, QOH-1 antibody, QOH1 antibody, QORL1_HUMAN antibody, Quinone oxidoreductase homolog 1 antibody, Quinone oxidoreductase like 1 antibody, Quinone oxidoreductase-like protein 1 antibody, Quinone reductase like 1 antibody, Zeta crystallin homolog antibody, Zeta-crystallin homolog antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	O95825
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Quinone oxidoreductase-like protein 1 protein (1-87AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:2000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western blot All lanes: CRYZL1 antibody at 12µg/ml
Lane 1: Hela whole cell lysate Lane 2: HepG2 whole cell lysate Lane 3: Jurkat whole cell lysate Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 39, 22 kDa Observed band size: 39 kDa