

IgG Antibody

PACO26501

Description

This IgG Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO26501
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	IgG
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Goat

Antibody Data

Isotype:	IgG
Uniprot:	-
Host Species:	Guinea pig
Purification:	>95%, Protein G purified
Immunogen:	Goat IgG heavy chain (Native Protein)
Immunogen Species:	Goat
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western blot All lanes: Goat IgG heavy chain antibody at 2µg/ml + Goat serum at 1: 100 Secondary Rabbit polyclonal to Guinea pig IgG at 1/15000 dilution
Predicted band size: 55 kDa Observed band size: 55 kDa