

Gstp1 Antibody

PACO27229

Description

This Gstp1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO27229
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Gstp1 antibody, GstpibGlutathione S-transferase P 1 antibody, Gst P1 antibody, EC 2.5.1.18 antibody, GST YF-YF antibody, GST class-pi antibody, GST-piB antibody, Preadipocyte growth factor antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Mouse, Human

Antibody Data

Isotype:	IgG
Uniprot:	P19157
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Mouse Glutathione S-transferase P 1 protein (2-210AA)
Immunogen Species:	Mus musculus (Mouse)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western blot All lanes: Gstp1 antibody at 2µg/ml
Lane 1: Mouse kidney tissue Lane 2: Jurkat whole cell lysate Lane 3: K562 whole cell lysate Lane 4: 293T whole cell lysate Lane 5: Hela whole cell lysate Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 24 kDa Observed band size: 24 kDa