

KRT38 Antibody

PACO27857

Description

This KRT38 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO27857
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	HA8 antibody, Hair keratin antibody, Hair keratin, type I Ha8 antibody, HHa8 antibody, HKA8 antibody, K38 antibody, Keratin 38 antibody, Keratin antibody, Keratin, hair, acidic, 8 antibody, Keratin, type I cuticular Ha8 antibody, Keratin-38 antibody, KRT38 antibody, KRT38_HUMAN antibody, KRTHA8 antibody, type I cuticular Ha8 antibody, type I Ha8 antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	O76015
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Keratin, type I cuticular Ha8 protein (2-456AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western blot All lanes: KRT38 antibody at 2µg/ml Lane 1: MCF-7 whole cell lysate Lane 2: Hela whole cell lysate Lane 3: HepG2 whole cell lysate Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 51 kDa Observed band size: 51 kDa