

MREG Antibody

PACO28330

Description

This MREG Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO28330
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	MREG antibody, DSU antibody, HDCGA21P antibody, Melanoregulin antibody, Dilute suppressor protein homolog antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse

Antibody Data

Isotype:	IgG
Uniprot:	Q8N565
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Melanoregulin protein (1-214AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western blot All lanes: MREG antibody at 8µg/ml + Mouse stomach tissue Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 25, 27 kDa Observed band size: 25, 33 kDa