

Tpsb2 Antibody

PACO30747

Description

This Tpsb2 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO30747
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Tpsb2 antibody, Mcpt6 antibody, Tryptase beta-2 antibody, Tryptase-2 antibody, EC 3.4.21.59 antibody, Mast cell protease 6 antibody, mMCP-6 antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	P21845
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Mouse Tryptase beta-2 protein (32-276AA)
Immunogen Species:	Mus musculus (Mouse)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:2000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western blot All lanes: Tpsb2 antibody at 2µg/ml Lane 1: Mouse heart tissue Lane 2: Rat heart tissue Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 31, 26 kDa Observed band size: 31 kDa