

ZNF300 Antibody

PACO31256

Description

This ZNF300 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO31256
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	DKFZp686B24204 antibody, Kruppel like zinc finger protein antibody, OTTHUMP00000160554 antibody, OTTHUMP00000206798 antibody, Zinc finger protein 300 antibody, ZN300_HUMAN antibody, ZNF300 antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q96RE9
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Zinc finger protein 300 protein (1-300AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-