

Adipoq Antibody

PACO32520

Description

This Adipoq Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO32520
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Adipoq antibody, Acdc antibody, Acrp30 antibody, Apm1 antibody, Adiponectin antibody, 30 kDa adipocyte complement-related protein antibody, Adipocyte complement-related 30 kDa protein antibody, ACRP30 antibody, Adipocyte antibody, C1q and collagen domain-containing protein antibody, Adipocyte-specific protein AdipoQ antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Mouse

Antibody Data

Isotype:	IgG
Uniprot:	Q60994
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Mouse Adiponectin protein (18-247AA)
Immunogen Species:	Mus musculus (Mouse)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-