

dgcQ Antibody

PACO32620

Description

This dgcQ Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO32620
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	dgcQ antibody, yedQ antibody, b1956 antibody, JW5832 antibody, Probable diguanylate cyclase DgcQ antibody, DGC antibody, EC 2.7.7.65 antibody, Cellulose synthesis regulatory protein antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Escherichia coli

Antibody Data

Isotype:	IgG
Uniprot:	P76330
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Escherichia coli Probable diguanylate cyclase DgcQ protein (381-564AA)
Immunogen Species:	Escherichia coli
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.


Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western Blot Positive WB detected in: Rosseta bacteria lysate at 80µg, 40µg, 20µg, 10µg All lanes: yedQ antibody at 2.5µg/ml Secondary Goat polyclonal to rabbit IgG at 1/50000 dilution predicted band size: 65 kDa observed band size: 65 kDa