

Ribosome-inactivating protein cucurmosin Antibody

PACO34042

Description

This Ribosome-inactivating protein cucurmosin Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO34042
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Ribosome-inactivating protein cucurmosin antibody, EC 3.2.2.22 antibody, rRNA N-glycosidase antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Cucurbita moschata

Antibody Data

Isotype:	IgG
Uniprot:	P84531
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Cucurbita moschata Ribosome-inactivating protein cucurmosin protein (1-244AA)
Immunogen Species:	Cucurbita moschata
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.


Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western Blot Positive WB detected in: Cucurbita moschata (40µg, 20µg, 10µg, 5µg) All lanes: rRNA N-glycosidase antibody at 4µg/ml Secondary Goat polyclonal to rabbit IgG at 1/50000 dilution Predicted band size: 28 kDa Observed band size: 28 kDa