

## Pka-C1 Antibody

PACO34146

### Description

---

This Pka-C1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	PACO34146
<b>Contents:</b>	50µg Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	-
<b>Synonyms:</b>	Pka-C1 antibody, CdkA antibody, DC0 antibody, CG4379 antibody, cAMP-dependent protein kinase catalytic subunit 1 antibody, PKA C antibody, EC 2.7.11.11 antibody, Protein kinase DC0 antibody
<b>Clone:</b>	Polyclonal
<b>Applications:</b>	<a href="#">ELISA</a>
<b>Conjugation:</b>	Non-conjugated
<b>Reactivity:</b>	Drosophila melanogaster

### Antibody Data

---

<b>Isotype:</b>	IgG
<b>Uniprot:</b>	P12370
<b>Host Species:</b>	Rabbit
<b>Purification:</b>	>95%, Protein G purified
<b>Immunogen:</b>	Recombinant Drosophila melanogaster cAMP-dependent protein kinase catalytic subunit 1 protein (2-353AA)
<b>Immunogen Species:</b>	Drosophila melanogaster
<b>Buffer:</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Form:</b>	Liquid

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

## Preparation & Storage

---

**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

---

Image	Description
-	-