

Mcpt4 Antibody

PACO34310

Description

This Mcpt4 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO34310
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Mast cell protease 4 antibody, Mcp 4 antibody, Mcpt4 antibody, MCPT4_MOUSE antibody, MMCP 4 antibody, MMCP 4A antibody, MMCP 4B antibody, MMCP-4 antibody, MMCP4 antibody, MMCP4A antibody, MMCP4B antibody, MSMCP antibody, Myonase antibody, Serosal mast cell protease antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Mouse

Antibody Data

Isotype:	IgG
Uniprot:	P21812
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Mouse Mast cell protease 4 protein (21-246AA)
Immunogen Species:	Mus musculus (Mouse)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-