

## Klk1b22 Antibody

PACO34314

### Description

---

This Klk1b22 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	PACO34314
<b>Contents:</b>	50µg Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	-
<b>Synonyms:</b>	Klk1b22 antibody, Klk-22 antibody, Klk22 antibody, Kallikrein 1-related peptidase b22 antibody, EC 3.4.21.35 antibody, Beta-NGF-endopeptidase antibody, Epidermal growth factor-binding protein type A antibody, EGF-BP A antibody, Glandular kallikrein K22 antibody, mGK-22 antibody, Nerve growth factor beta chain endopeptidase antibody, Tissue kallikrein 22 antibody
<b>Clone:</b>	Polyclonal
<b>Applications:</b>	<b>ELISA</b>
<b>Conjugation:</b>	Non-conjugated
<b>Reactivity:</b>	Mouse

### Antibody Data

---

<b>Isotype:</b>	IgG
<b>Uniprot:</b>	P15948
<b>Host Species:</b>	Rabbit
<b>Purification:</b>	>95%, Protein G purified
<b>Immunogen:</b>	Recombinant Mouse Kallikrein 1-related peptidase b22 protein (25-259AA)
<b>Immunogen Species:</b>	Mus musculus (Mouse)
<b>Buffer:</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Form:** Liquid

## Preparation & Storage

---

**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

<b>Recommended Dilutions:</b>	<b>Application</b>	<b>Recommended Dilution</b>
	-	-

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

---

<b>Image</b>	<b>Description</b>
-	-