

Reg3b Antibody

PACO34350

Description

This Reg3b Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO34350
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Reg3b antibody, Pap antibody, Pap1 antibody, Reg2 antibody, Regenerating islet-derived protein 3-beta antibody, REG-3-beta antibody, Pancreatitis-associated protein 1 antibody, Peptide 23 antibody, REG-2 antibody, Regenerating islet-derived protein III-beta antibody, Reg III-beta) [Cleaved into: Regenerating islet-derived protein 3-beta 16.5 kDa form, Regenerating islet-derived protein 3-beta 15 kDa form] antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Rat

Antibody Data

Isotype:	IgG
Uniprot:	P25031
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Rat Regenerating islet-derived protein 3-beta protein (27-175AA)
Immunogen Species:	Rattus norvegicus (Rat)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-