

## Pectate lyase 5 Antibody

PACO34386

### Description

---

This Pectate lyase 5 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	PACO34386
<b>Contents:</b>	50µg Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	-
<b>Synonyms:</b>	Pectate lyase 5 antibody, EC 4.2.2.2 antibody, Antigen Amb a I antibody, Antigen E antibody, AgE antibody, Pollen allergen Amb a 1.1 antibody, allergen Amb a 1.1 antibody
<b>Clone:</b>	Polyclonal
<b>Applications:</b>	<b>ELISA</b>
<b>Conjugation:</b>	Non-conjugated
<b>Reactivity:</b>	Ambrosia artemisiifolia

### Antibody Data

---

<b>Isotype:</b>	IgG
<b>Uniprot:</b>	P27759
<b>Host Species:</b>	Rabbit
<b>Purification:</b>	>95%, Protein G purified
<b>Immunogen:</b>	Recombinant Ambrosia artemisiifolia Pectate lyase 5 protein (26-396AA)
<b>Immunogen Species:</b>	Ambrosia artemisiifolia (Short ragweed)
<b>Buffer:</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Form:</b>	Liquid

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

## Preparation & Storage

---

**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

---

Image	Description
-	-