

tetX Antibody

PACO34670

Description

This tetX Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO34670
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	tetX antibody, CTC_p60 antibody, Tetanus toxin antibody, EC 3.4.24.68 antibody, Tentoxylisin) [Cleaved into: Tetanus toxin light chain antibody, Tetanus toxin chain L), Tetanus toxin heavy chain antibody, Tetanus toxin chain H)] antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Clostridium tetani

Antibody Data

Isotype:	IgG
Uniprot:	P04958
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Clostridium tetani Tetanus toxin protein (123-573AA)
Immunogen Species:	Clostridium tetani (strain Massachusetts / E88)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-