

surA Antibody

PACO34790

Description

This surA Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO34790
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	surA antibody, b0053 antibody, JW0052 antibody, Chaperone SurA antibody, Peptidyl-prolyl cis-trans isomerase SurA antibody, PPlase SurA antibody, EC 5.2.1.8 antibody, Rotamase SurA antibody, Survival protein A antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Escherichia coli

Antibody Data

Isotype:	IgG
Uniprot:	P0ABZ6
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Escherichia coli Chaperone SurA protein (21-428AA)
Immunogen Species:	Escherichia coli
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western blot All lanes: surA antibody at 12µg/ml Lane 1: 293T whole cell lysate Lane 2: DH5α whole cell lysate Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 48 kDa Observed band size: 48 kDa