

crp Antibody

PACO34794

Description

This crp Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO34794
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	crp antibody, cap antibody, csm antibody, b3357 antibody, JW5702cAMP-activated global transcriptional regulator CRP antibody, Catabolite activator protein antibody, CAP antibody, Catabolite gene activator antibody, cAMP receptor protein antibody, CRP antibody, cAMP regulatory protein antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Escherichia coli

Antibody Data

Isotype:	IgG
Uniprot:	P0ACJ8
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Escherichia coli cAMP-activated global transcriptional regulator CRP protein (1-210AA)
Immunogen Species:	Escherichia coli (strain K12)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-