

glaB1 Antibody

PACO35018

Description

This glaB1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO35018
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	glaB1 antibody, BVU_0104 antibody, Alpha-1,3-galactosidase B antibody, EC 3.2.1.n1 antibody, EC 3.2.1.n2 antibody, Exo-alpha-galactosidase B1 antibody, EC 3.2.1.22 antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Bacteroides vulgatus

Antibody Data

Isotype:	IgG
Uniprot:	A6KWM0
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Bacteroides vulgatus Alpha-1,3-galactosidase B protein (19-592AA)
Immunogen Species:	Bacteroides vulgatus (strain ATCC 8482 / DSM 1447 / JCM 5826 / NBRC 14291 / NCTC 11154)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-