

## MRJP1 Antibody

PACO35110

### Description

---

This MRJP1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	PACO35110
<b>Contents:</b>	50µg Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	-
<b>Synonyms:</b>	MRJP1 antibody, Major royal jelly protein 1 antibody, MRJP-1 antibody, MRJP1 antibody, 56-kDa protein 4 antibody, p56kP-4 antibody, Apalbumin 1 antibody, Apisin subunit MRJP1 antibody, Bee-milk protein antibody, Royalactin) [Cleaved into: Jellein-1 antibody, Jelleine-I), Jellein-2 antibody, Jelleine-II), Jellein-4 antibody, Jelleine-IV)] antibody
<b>Clone:</b>	Polyclonal
<b>Applications:</b>	<b>ELISA</b>
<b>Conjugation:</b>	Non-conjugated
<b>Reactivity:</b>	Apis mellifera

### Antibody Data

---

<b>Isotype:</b>	IgG
<b>Uniprot:</b>	O18330
<b>Host Species:</b>	Rabbit
<b>Purification:</b>	>95%, Protein G purified
<b>Immunogen:</b>	Recombinant Apis mellifera Major royal jelly protein 1 protein (20-432AA)
<b>Immunogen Species:</b>	Apis mellifera (Honeybee)
<b>Buffer:</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Form:</b>	liquid

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

## Preparation & Storage

---

**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

---

Image	Description
-	-