

KCTD19 Antibody

PACO35814

Description

This KCTD19 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO35814
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	BTB/POZ domain-containing protein KCTD19 antibody, FLJ40162 antibody, KCD19_HUMAN antibody, KCTD19 antibody, potassium channel tetramerisation domain containing 19 antibody, Testicular tissue protein Li 101 antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q17RG1
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human BTB/POZ domain-containing protein KCTD19 protein (490-751AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human liver cancer using PACO35814 at dilution of 1:100

Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO35814 at dilution of 1:100