

ANKDD1A Antibody

PACO36002

Description

This ANKDD1A Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO36002
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	AKD1A_HUMAN antibody, ANKDD1A antibody, Ankyrin repeat and death domain containing 1A antibody, Ankyrin repeat and death domain-containing protein 1A antibody, FLJ18078 antibody, FLJ25870 antibody, MGC120305 antibody, MGC120306 antibody, MGC120307 antibody, OTTHUMP00000164083 antibody, OTTHUMP00000164131 antibody, OTTHUMP00000220561 antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q495B1
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Ankyrin repeat and death domain-containing protein 1A protein (272-490AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol


Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human testis tissue using PACO36002 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO36002 at dilution of 1:100