

C8orf34 Antibody

PACO36030

Description

This C8orf34 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO36030
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	C8orf34Uncharacterized protein C8orf34 antibody, Protein VEST-1 antibody
Clone:	Polyclonal
Applications:	ELISA WB IHC
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse

Antibody Data

Isotype:	IgG
Uniprot:	Q49A92
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Uncharacterized protein C8orf34 protein (1-300AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:5000
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol


Validation Data

Image




Description

Western Blot All lanes: C8orf34 antibody at 5µg/ml + Mouse skeletal muscle tissue Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 60, 48, 33, 23 kDa Observed band size: 60 kDa



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO36030 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skin tissue using PACO36030 at dilution of 1:100