

fabF Antibody

PACO36202

Description

This fabF Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO36202
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	fabF antibody, SACOL09883-oxoacyl-[acyl-carrier-protein] synthase 2 antibody, EC 2.3.1.179 antibody, 3-oxoacyl-[acyl-carrier-protein] synthase II antibody, Beta-ketoacyl-ACP synthase II antibody, KAS II antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Staphylococcus aureus

Antibody Data

Isotype:	IgG
Uniprot:	Q5HHA1
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Staphylococcus aureus 3-oxoacyl-[acyl-carrier-protein] synthase 2 protein (1-414AA)
Immunogen Species:	Staphylococcus aureus (strain COL)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-