

KRT79 Antibody

PACO36426

Description

This KRT79 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO36426
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	KRT79 antibody, K6L antibody, KB38 antibody, KRT6L antibody, Keratin antibody, type II cytoskeletal 79 antibody, Cytokeratin-79 antibody, CK-79 antibody, Keratin-6-like antibody, Keratin-6L antibody, Keratin-79 antibody, K79 antibody, Type-II keratin Kb38 antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse

Antibody Data

Isotype:	IgG
Uniprot:	Q5XKE5
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Keratin, type II cytoskeletal 79 protein (1-300AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western blot All lanes: KRT79 antibody at 1µg/ml Lane 1: Mouse lung tissue Lane 2: A549 whole cell lysate Lane 3: 293T whole cell lysate Lane 4: Mouse gonadal tissue Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 58 kDa Observed band size: 58 kDa