

LGALS9C Antibody

PACO36526

Description

This LGALS9C Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO36526
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	LGALS9C antibody, Galectin-9C antibody, Gal-9C antibody, Galectin-9-like protein B antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse

Antibody Data

Isotype:	IgG
Uniprot:	Q6DKI2
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Galectin-9C protein (1-356AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western blot All lanes: LGALS9C antibody at 4µg/ml
Lane 1: Mouse stomach tissue Lane 2: Mouse liver tissue
Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 40 kDa Observed band size: 22 kDa