

MRNIP Antibody

PACO36826

Description

This MRNIP Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO36826
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	MRNIP antibody, C5orf45 antibody, MRN complex-interacting protein antibody, MRN-interacting protein antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse

Antibody Data

Isotype:	IgG
Uniprot:	Q6NTE8
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human MRN complex-interacting protein (1-200AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.


Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western blot All lanes: MRNIP antibody at 4µg/ml Lane 1: Mouse kidney tissue Lane 2: Mouse heart tissue Lane 3: Hela whole cell lysate Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 38, 30, 32 kDa Observed band size: 38 kDa